

Intrauterin Insemination

CATHETERS (IUI)

www.kitazato.co.jp



Intrauterin Insemination CATHETERS (IUI)

- 4 types of rigidities: Hard, Ultimate, Intermediate and Soft
- Stylet cannula specially designed for single step IUI
- Two side holes for better sperm dispersion
- Smooth and rounded tip

IUI Catheters



Intrauterine Insemination Catheters (IUI) - Long Type

| Order Number | Catheter (Fr) | Catheter Length (cm) | Total volume minus hub | Rigidity | Hub Color | Quantity |
|--------------|---------------|----------------------|------------------------|--------------|------------|----------|
| 131181 | 6 | 18 | 0.03mL | Hard | Dark gray | 50 pcs |
| 131186 | 6 | 18 | 0.03mL | Ultimate | Gray | 50 pcs |
| 131182 | 6 | 18 | 0.10mL | Intermediate | Light gray | 50 pcs |
| 131184 | 5.2 | 18 | 0.10mL | Soft | White | 50 pcs |

Stylet Cannula for Soft and Intermediate type.

Order Number Total volume minus hub used with Soft and Intermediate type (131182 & 131184)

| | |
|--------|--------|
| 140198 | 0.03mL |
|--------|--------|

Intrauterine Insemination Catheters (IUI) - Short Type

| Order Number | Catheter (Fr) | Catheter Length (cm) | Total volume minus hub | Rigidity | Hub Color | Quantity |
|--------------|---------------|----------------------|------------------------|--------------|------------|----------|
| 131101 | 6 | 10 | 0.015mL | Hard | Dark gray | 50 pcs |
| 131106 | 6 | 10 | 0.015mL | Ultimate | Gray | 50 pcs |
| 131102 | 6 | 10 | 0.05mL | Intermediate | Light gray | 50 pcs |
| 131104 | 5.2 | 10 | 0.05mL | Soft | White | 50 pcs |

Stylet Cannula for Soft and Intermediate type.

Order Number Total volume minus hub used with Soft and Intermediate type (131102 & 131104)

| | |
|--------|---------|
| 140118 | 0.015mL |
|--------|---------|

Composition

- Polyvinylchloride
- Stainless Steel SUS304

Quality Control

- MEA (Mouse Embryo Assay): One-cell assay $\geq 80\%$ after 96 hours
- Endotoxin < 20 EU/device
- SAL 10^{-6}
- Sterility Test (Bacteria, Fungi)
- HSSA (Human Sperm Survival Assay): $\geq 70\%$ motility at 24 hours
- Cytotoxicity Test
- Intracutaneous Reactivity Test
- Sensitization Test



IUI **Soft / Intermediate**



IUI **Hard** with Built-in Stylet



Dibimed

Manufactured by:
Kitazato Medical Co. Ltd.
1-1-8 Shibadaimon, Minato-ku
Tokyo 105-0012 JAPAN
Tel +81-3-3434-2731
Fax +81-3-3434-2732
trading@kitazato.co.jp

Distributed by:
Biomedical Supply, S.L.
Luis Buñuel, 1 - Of. Pta. 4
46015 Valencia SPAIN
Tel (+34) 963 056 395
Fax (+34) 963 056 396
info@dibimed.com