

Embryo Transfer

CATHETERS



www.kitazato.co.jp



Supported Embryo Transfer CATHETERS

- Echo marked tip
- 2 choices of end: straight or bulb type
- Obturators and Stylets available
- 30° precurved catheters

The range of Kitazato supported catheters includes inner catheters in 3 and 4 French sizes. We have also developed the extra-long 40cm catheter for the transvaginal ultrasound guided embryo transfer.



Composition (Supported ET Catheters)

Catheter: Silicone, Stainless Steel SUS304
Guide: 12 Nylon

Supported Embryo Transfer Straight Type Catheters



Catheter and Guide



Stylet



Catheter and Guide with Obturator

Supported Embryo Transfer Catheters - Straight Type

Order Number	Catheter Length (cm)	Box	Catheter + Guide	Obturator	Stylet
214322	22	10 pcs	x	--	--
214325	25	10 pcs	x	--	--
224322	22	10 pcs	x	x	--
224325	25	10 pcs	x	x	--
244017	--	10 pcs	--	--	x
244020	--	10 pcs	--	--	x

Routine Embryo Transfer CATHETERS

- Soft Catheter to minimize trauma
- Keep memory when giving shape
- Rounded tip for smoother insertion
- Easy hub to introduce the catheter

Routine Embryo Transfer Catheters						
Order Number	Catheter Length (cm)	Box	Catheter + Guide	Obturator	Stylet	
215018	18	10 pcs	x	--	--	
215023	23	10 pcs	x	--	--	


Kitazato ET Catheters are designed under the highest standards of quality and the suggestions of the most recognized gynecologists in order to maximize the embryo implantation.

Composition

(Routine ET Catheters)
 Catheter: Polyurethane
 Guide: Fluorocarbon resin

Quality Control

- MEA (Mouse Embryo Assay): One-cell assay \geq 80% after 96 hours
- Endotoxin: <20EU/device
- SAL10⁻⁶
- Sterility Test (Bacteria, Fungi)
- Cytotoxicity Test
- Intracutaneous Reactivity Test
- Sensitization Test



Supported Embryo Transfer Bulb Type Catheters



Catheter and Guide



Stylet



Catheter and Guide with Obturator

Supported Embryo Transfer Catheters - Bulb Type

Order Number	Catheter Length (cm)	Box	Catheter + Guide	Obturator	Stylet
213322	22	10 pcs	x	--	--
213325	25	10 pcs	x	--	--
223322	22	10 pcs	x	x	--
223325	25	10 pcs	x	x	--
243017	--	10 pcs	--	--	x
243020	--	10 pcs	--	--	x

Trans Vaginal Ultrasound Guided Embryo Transfer Catheters



Catheter and Guide with Obturator (Stylet available)

Supported Embryo Transfer Catheters - Bulb Type

Order Number	Catheter Length (cm)	Box	Catheter + Guide	Obturator	Stylet
223340	40	10 pcs	x	x	--
243020	--	10 pcs	--	--	x

Clinical Reference

Bodri, D. - Transvaginal versus transabdominal ultrasound guidance for embryo transfer in donor oocyte recipients: a randomized clinical trial, Fertility & Sterility 2011.



Dibimed

Manufactured by:
Kitazato Medical Co. Ltd.
2-5-12 Hongo, Bunkyo-Ku
Tokyo 113-0033 JAPAN
Tel +81-3-5800-6055
Fax 81-3-5800-6056
trading@kitazato.co.jp

Distributed by:
Biomedical Supply, S.L.
Luis Buñuel, 1 - Of. Pta. 4
46015 Valencia SPAIN
Tel (+34) 963 056 395
Fax (+34) 963 056 396
info@dibimed.com