

Oocyte Retrieval **NEEDLES**

www.kitazato.co.jp



Oocyte Retrieval NEEDLES

Single Lumen Oocyte Aspiration

- Echo marked tip for better positioning
- Triple cut tip to minimize trauma
- Best grip design for firm control
- Silicon bung for safety in test tubes
- Thin wall needles for smoother aspiration

Kitazato Needles have been designed under the highest standards of quality in order to maximize the quality of oocytes retrieved in the minimum operating time and patient pain.

The range of Kitazato is the widest of all IVF manufacturers, from 16G to 21G. The reduction of the needle tip allows less pain, less tissue damage and less trauma without affecting the number of oocytes collected, their quality or the clinical pregnancy rate.




Composition

Stainless Steel SUS304

Quality Control

- MEA (Mouse Embryo Assay): One-cell assay \geq 80% after 96 hours
- Endotoxin <20 EU/device
- SAL 10^{-6}
- Sterility Test (Bacteria, Fungi)
- Cytotoxicity Test
- Intracutaneous Reactivity Test
- Sensitization Test

Aspiration Sets for Pump



Triple cut needle tip
Echo marked tip



Ergonomic grip



Aspiration Sets for Pump

Order Number	Gauge	Length (cm)	Aspiration Line (cm)*	Box
326350	16	35	80	10 pcs
327350	17	35	80	10 pcs
328350	18	35	80	10 pcs
329350	19	35	80	10 pcs
320300	20	30	80	10 pcs
320350	20	35	80	10 pcs
321300	21	30	80	10 pcs

*Other lengths of aspiration line available

Oocyte Retrieval NEEDLES

Aspiration Sets for Syringe



Triple cut needle tip
Echo marked tip



Ergonomic grip



Aspiration Sets for Syringe

Order Number	Gauge	Length (cm)	Aspiration Line (cm)*	Box
367350	17	35	60	10 pcs

*Other lengths of aspiration line available

Aspiration Sets with Luer Lock



Luer Lock for aspiration line/syringe




Aspiration Sets with Luer Lock

Order Number	Gauge	Length (cm)	Aspiration Line (cm)	Box
316001	16	32.5	--	10 pcs
317001	17	32.5	--	10 pcs
318001	18	32.5	--	10 pcs
319001	19	32.5	--	10 pcs
310001	20	32.5	--	10 pcs

Oocyte Retrieval NEEDLES

Aspiration Sets with Stopcock



3 way stopcock



Silicon bung



Aspiration Sets with Stopcock

Order Number	Gauge	Length (cm)	Aspiration Line (cm)*	Box
357353	17	35	80	20 pcs

*Other lengths of aspiration line available

Aspiration Set for IVM



20G aspiration needle inside a 17G needle



Aspiration Set for IVM

Order Number	Gauge	Length (cm)	Aspiration Line (cm)	Box
347325	17	27.5	30	10 pcs
	20	35		



Dibimed

Manufactured by:
Kitazato Medical Co. Ltd.
2-5-12 Hongo, Bunkyo-Ku
Tokyo 113-0033 JAPAN
Tel +81-3-5800-6055
Fax 81-3-5800-6056
trading@kitazato.co.jp

Distributed by:
Biomedical Supply, S.L.
Luis Buñuel, 1 - Of. Pta. 4
46015 Valencia SPAIN
Tel (+34) 963 056 395
Fax (+34) 963 056 396
info@dibimed.com