

SUPPORTED AND
FULLY ECHOGENIC
ET CATHETER

RYDEN™

Experience the difference.
Clarity and control every time.



New

- **Fully echogenic** catheter with echo-tip on the guide
- **Outer metallic support** for more control
- **Bulb tip** for atraumatic transfer

Ultrasound image provided by Atsumi Yoshida, MD, PhD,
Kiba Park Clinic, JAPAN

KITAZATO®

PRODUCT REFERENCES

RYDEN™ ET Catheter


| Order Number | Description | Diameter (Fr) | Length (cm) |
|--------------|---|---------------|-------------|
| 213312 | Echogenic Supported Catheter with Guide | 3 | 22 |
| 213355 | Echogenic Supported Catheter with Guide | 3 | 25 |
| 223312 | Echogenic Supported Catheter and Guide with Obturator | 3 | 22 |
| 223355 | Echogenic Supported Catheter and Guide with Obturator | 3 | 25 |

| Order Number | Description |
|--------------|-------------------------------------|
| 244017 | Stylet for guide with 22cm catheter |
| 244020 | Stylet for guide with 25cm catheter |



Catheter with metallic support that generates rigidity for a better control of the catheter

Echogenicity along the whole catheter, which will provide better visualization along the way to the uterine cavity



Bulb-tip with Echogenic ring that enhances a safe, smooth and atraumatic pass-way through the cervix

Pre-curved guide with memory to give an easy pass along



Large and user-friendly opening in the handle for insertion of the preloaded catheter into the guide